



HeraCeram® Zirkonia 750

Ceramic for veneering zirconia and lithium disilicate frameworks



● **Optically perfect – with high-purity quartz glass**

● **Technically perfect – with a stabilised leucite structure (SLS)**

● **Simply perfect – with a consistent processing philosophy**



IMPORTANT: The firing temperatures quoted here are guidelines only. Deviations may occur due to differences in furnace performance and may have to be compensated for if necessary.

General Firing programmes

For Zirconiumdioxide-Frameworks							For Lithium-Disilicate-Frameworks						
Adhesive Margin HM Dentine1 Dentine2 Glaze Correction COR/LM							Adhesive Margin HM Dentine1 Dentine2 Glaze Correction COR/LM						
Preheating or staring temperature[°C]	500	500	500	500	500	500	400	400	400	400	400	400	
Predrying and pre- heating time [min]	6	6	6	6	5	5	6	6	6	6	5	5	
Heat rate [°C/min]	75	60	60	60	60	60	50	50	50	50	60	50	
Final temp. [°C]	800	780	750	745	725	720	750	780	750	745	725	715	
Holding time [°C]	2	1	1	1	0.5-1	1	1	1	1	1	0.5-1	1	
Vacuum start [°C]	500	500	500	500	–	500	400	400	400	400	–	400	
Vacuum tim stop [°C]	800	780	750	745	–	720	750	780	750	745	–	715	

Everyday A – D Shades

Reliable reproduction
of classic A – D
shades

●Dentine Incisal set and/or ●Chroma dentine set

Optional:

●Stains set universal

Personalised

Patient specific
shade adjustment

●Dentine Incisal set and/or ●Chroma dentine set

Optional:

●Stains set universal

● Increaser set
● Enhancer set



Our starter set: First-Touch set (A2; A3; A3,5)

Matrix

Truly natural for highly
individualised restorations.

●Dentine Incisal set and/or ●Chroma dentine set

Optional:

●Stains set universal

● Matrix set



Special

Special applications

Margin set

Bleach
shades

Gingiva set